

Insulin-like growth factor-binding protein 5

Catalog No: #AP89409

Package Size: #AP89409-1 50ug #AP89409-2 100ug #AP89409-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor-binding protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IBP-5
Accession No.	Uniprot ID:Q05717
Uniprot	Q05717
GeneID	404185;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IGFBP5

Note: This product is for in vitro research use only