

OR5T1 Antibody

Catalog No: #44614

Package Size: #44614-1 50ul #44614-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OR5T1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	OR5T1 Antibody detects endogenous levels of total OR5T1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR5T1
Target Name	OR5T1
Other Names	Olfactory receptor 5T1 Olfactory receptor OR11-179,OR5T1P
Accession No.	Swiss-Prot#: Q8NG75NCBI Gene ID:390155
Uniprot	Q8NG75
GeneID	390155;
Calculated MW	37 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis OR5T1 using 293 whole cell lysates

Note: This product is for in vitro research use only