

Lysophospholipid acyltransferase 5

Catalog No: #AP89427



Package Size: #AP89427-1 50ug #AP89427-2 100ug #AP89427-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Lysophospholipid acyltransferase 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	LPLAT 5,1-acylglycerophosphocholine O-acyltransferase,Lysophosphatidylcholine acyltransferase 3,Lyso-PC acyltransferase 3,Membrane-bound O-acyltransferase domain-containing protein 5,O-acyltransferase domain-containing protein 5,MBOAT5,OACT5
Accession No.	Uniprot ID:Q3SZL3
Uniprot	Q3SZL3
GeneID	515361;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LPCAT3

Note: This product is for in vitro research use only