

Macrophage migration inhibitory factor

Catalog No: #AP89431

Package Size: #AP89431-1 50ug #AP89431-2 100ug #AP89431-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Macrophage migration inhibitory factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MIF,Glutathione-binding 13 kDa protein,L-dopachrome isomerase,L-dopachrome tautomerase,Phenylpyruvate tautomerase
Accession No.	Uniprot ID:P30904
Uniprot	P30904
GeneID	103694877;81683;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mif

Note: This product is for in vitro research use only