

Metalloproteinase inhibitor 3

Catalog No: #AP89494



Package Size: #AP89494-1 50ug #AP89494-2 100ug #AP89494-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Metalloproteinase inhibitor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tissue inhibitor of metalloproteinases 3,TIMP-3,Sun,Timp-3
Accession No.	Uniprot ID:P39876
Uniprot	P39876
GeneID	21859;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Timp3

Note: This product is for in vitro research use only