

Vitamin D3 receptor

Catalog No: #AP89503

Package Size: #AP89503-1 50ug #AP89503-2 100ug #AP89503-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vitamin D3 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	VDR,1,25-dihydroxyvitamin D3 receptor,Nuclear receptor subfamily 1 group I member 1,Nr1i1
Accession No.	Uniprot ID:P48281
Uniprot	P48281
GeneID	22337;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Vdr

Note: This product is for in vitro research use only