

Rat Bmp4/Bone morphogenetic protein 4 ELISA Kit

Catalog No: #EK19176

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bmp4/Bone morphogenetic protein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bone morphogenetic protein 4
Target Name	Bone morphogenetic protein 4
Other Names	BMP-4, Bone morphogenetic protein 2B, BMP-2B, Bmp-4, Dvr-4
Accession No.	Q06826
Uniprot	Q06826

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmp4

Note: This product is for in vitro research use only