

## Chicken CNTF/Ciliary neurotrophic factor ELISA Kit

Catalog No: #EK19185

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken CNTF/Ciliary neurotrophic factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Ciliary neurotrophic factor
Target Name	Ciliary neurotrophic factor
Other Names	CNTF, Growth-promoting activity, GPA, CNF, GPA
Accession No.	Q02011
Uniprot	Q02011
GeneID	396488;

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000 pg/mL Sensitivity: 15.6 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene\_full\_name: CNTF

Note: This product is for in vitro research use only