



# Mouse Csf2/Granulocyte-macrophage colony-stimulating factor ELISA Kit

Catalog No: #EK19196

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Csf2/Granulocyte-macrophage colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Granulocyte-macrophage colony-stimulating factor
Target Name	Granulocyte-macrophage colony-stimulating factor
Other Names	GM-CSF, Colony-stimulating factor, CSF, Csfgm
Accession No.	P01587
Uniprot	P01587
GeneID	12981;

## Application Details

Assay Type: SandwichDetect Range:1.56-100pg/mL Sensitivity:0.66pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Csf2

Note: This product is for in vitro research use only