

Bovine ANGPT4/Angiopoietin-4 ELISA Kit

Catalog No: #EK19346



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Bovine ANGPT4/Angiopoietin-4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Angiopoietin-4
Target Name	Angiopoietin-4
Other Names	ANG-4
Accession No.	Q24K15
Uniprot	Q24K15
GeneID	617915;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.157ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:ANGPT4

Note: This product is for in vitro research use only