

# Human CELA2A/Chymotrypsin-like elastase family member 2A ELISA Kit

Catalog No: #EK19453



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CELA2A/Chymotrypsin-like elastase family member 2A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Chymotrypsin-like elastase family member 2A
Target Name	Chymotrypsin-like elastase family member 2A
Other Names	Elastase-2A,ELA2A
Accession No.	P08217
Uniprot	P08217
GeneID	63036;

## Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.403ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:CELA2A

Note: This product is for in vitro research use only