

Bovine ARG1/Arginase-1 ELISA Kit

Catalog No: #EK19476

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine ARG1/Arginase-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Arginase-1
Target Name	Arginase-1
Other Names	Liver-type arginase, Type I arginase
Accession No.	Q2KJ64
Uniprot	Q2KJ64
GeneID	513608;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 41.7pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ARG1

Note: This product is for in vitro research use only