

Rat Abcb9/ATP-binding cassette sub-family B member 9 ELISA Kit

Catalog No: #EK19511



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Abcb9/ATP-binding cassette sub-family B member 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat ATP-binding cassette sub-family B member 9
Target Name	ATP-binding cassette sub-family B member 9
Other Names	ATP-binding cassette transporter 9,ABC transporter 9 protein,TAP-like protein,TAPL
Accession No.	Q9QYJ4
Uniprot	Q9QYJ4
GeneID	63886;

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.4ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Abcb9

Note: This product is for in vitro research use only