



Mouse Cndp2/Cytosolic non-specific dipeptidase ELISA Kit

Catalog No: #EK19513

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cndp2/Cytosolic non-specific dipeptidase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cytosolic non-specific dipeptidase
Target Name	Cytosolic non-specific dipeptidase
Other Names	CNDP dipeptidase 2,Glutamate carboxypeptidase-like protein 1,Cn2
Accession No.	Q9D1A2
Uniprot	Q9D1A2
GeneID	66054;

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.401ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cndp2

Note: This product is for in vitro research use only