

Rat Cast/Calpastatin ELISA Kit

Catalog No: #EK19526



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Cast/Calpastatin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Calpastatin
Target Name	Calpastatin
Other Names	Calpain inhibitor
Accession No.	P27321
Uniprot	P27321
GeneID	25403;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:35pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cast

Note: This product is for in vitro research use only