

Pig CCL25/C-C motif chemokine 25 ELISA Kit

Catalog No: #EK19548



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Pig CCL25/C-C motif chemokine 25 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig C-C motif chemokine 25
Target Name	C-C motif chemokine 25
Other Names	Chemokine ligand 25,Small-inducible cytokine A25
Accession No.	Q4PR21
Uniprot	Q4PR21
GeneID	448799;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CCL25

Note: This product is for in vitro research use only