

# Rat Aif1/Allograft inflammatory factor 1 ELISA Kit

Catalog No: #EK19551



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                            |
|--------------------|--------------------------------------------------------------------------------------------|
| Product Name       | Rat Aif1/Allograft inflammatory factor 1 ELISA Kit                                         |
| Brief Description  | ELISA Kits                                                                                 |
| Applications       | ELISA                                                                                      |
| Species Reactivity | Rt                                                                                         |
| Specificity        | Natural and recombinant Rat Allograft inflammatory factor 1                                |
| Target Name        | Allograft inflammatory factor 1                                                            |
| Other Names        | AIF-1,Ionized calcium-binding adapter molecule 1,MRF-1,Microglia response factor,Iba1,Mrf1 |
| Accession No.      | P55009                                                                                     |
| Uniprot            | P55009                                                                                     |
| GeneID             | 29427;                                                                                     |

## Application Details

Assay Type: SandwichDetect Range:15.6-1000pg/mL Sensitivity:7.82pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Aif1

Note: This product is for in vitro research use only