

Mouse Ar/Androgen receptor ELISA Kit

Catalog No: #EK19552

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Ar/Androgen receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Androgen receptor
Target Name	Androgen receptor
Other Names	Dihydrotestosterone receptor,Nuclear receptor subfamily 3 group C member 4,Nr3c4
Accession No.	P19091
Uniprot	P19091
GeneID	11835;

Application Details

Assay Type: SandwichDetect Range:31.2-2000pg/mL Sensitivity:12pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Ar

Note: This product is for in vitro research use only