

Rat Ahcy/Adenosylhomocysteinase ELISA Kit

Catalog No: #EK19564

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ahcy/Adenosylhomocysteinase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Adenosylhomocysteinase
Target Name	Adenosylhomocysteinase
Other Names	AdoHcyase, S-adenosyl-L-homocysteine hydrolase
Accession No.	P10760
Uniprot	P10760
GeneID	29443;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.68pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ahcy

Note: This product is for in vitro research use only