

## Mouse Cdh5/Cadherin-5 ELISA Kit

Catalog No: #EK19606

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cdh5/Cadherin-5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cadherin-5
Target Name	Cadherin-5
Other Names	Vascular endothelial cadherin, VE-cadherin, CD144
Accession No.	P55284
Uniprot	P55284
GeneID	12562;

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.33ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cdh5

Note: This product is for in vitro research use only