

Rat Cib1/Calcium and integrin-binding protein 1 ELISA Kit

Catalog No: #EK19612

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Cib1/Calcium and integrin-binding protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Calcium and integrin-binding protein 1
Target Name	Calcium and integrin-binding protein 1
Other Names	CIB,Calmyrin,DNA-PKcs-interacting protein,Kinase-interacting protein,KIP,Cib,Kip,Prkdcip
Accession No.	Q9R010
Uniprot	Q9R010
GeneID	81823;

Application Details

Assay Type: SandwichDetect Range:0.625-40ng/mL Sensitivity:0.267 ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cib1

Note: This product is for in vitro research use only