

Rat Chat/Choline O-acetyltransferase ELISA Kit

Catalog No: #EK19644



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Rat Chat/Choline O-acetyltransferase ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat Choline O-acetyltransferase |
| Target Name | Choline O-acetyltransferase |
| Other Names | CHOACTase |
| Accession No. | P32738 |
| Uniprot | P32738 |

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.438ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Chat

Note: This product is for in vitro research use only