

Mouse Casq1/Calsequestrin-1 ELISA Kit

Catalog No: #EK19695



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Casq1/Calsequestrin-1 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Calsequestrin-1 |
| Target Name | Calsequestrin-1 |
| Other Names | Calmitine,Calsequestrin, skeletal muscle isoform |
| Accession No. | O09165 |
| Uniprot | O09165 |
| GeneID | 12372; |

Application Details

Assay Type: SandwichDetect Range:0.078-5ng/mL Sensitivity:0.03ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Casq1

Note: This product is for in vitro research use only