

# Mouse Adcyap1r1/Pituitary adenylate cyclase-activating polypeptide type I receptor ELISA Kit

Catalog No: #EK19709

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Adcyap1r1/Pituitary adenylate cyclase-activating polypeptide type I receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Pituitary adenylate cyclase-activating polypeptide type I receptor
Target Name	Pituitary adenylate cyclase-activating polypeptide type I receptor
Other Names	PACAP type I receptor
Accession No.	P70205
Uniprot	P70205
GeneID	11517;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 39.7 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene\_full\_name: Adcyap1r1

Note: This product is for in vitro research use only