

UKHC antibody

Catalog No: #23030

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	UKHC antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 901 and 963 of Human KIF5B
Target Name	UKHC
Accession No.	Swiss-Prot:P33176Gene ID:3799
Uniprot	P33176
GeneID	3799;
Concentration	0.4mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 110kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Images



Sample (30 ug of whole cell lysate)


A: 293T

B: A431

C: H1299

15% SDS PAGE

Primary antibody diluted at 1: 1500



Immunohistochemical analysis of paraffin-embedded Breast ca, using UKHC antibody at 1: 500 dilution.

Note: This product is for in vitro research use only