



Mouse Angptl3/Angiopoietin-related protein 3 ELISA Kit

Catalog No: #EK19792

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Mouse Angptl3/Angiopoietin-related protein 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiopoietin-related protein 3
Target Name	Angiopoietin-related protein 3
Other Names	Angiopoietin-like protein 3
Accession No.	Q9R182
Uniprot	Q9R182
GeneID	30924;

Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.414ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Angptl3

Note: This product is for in vitro research use only