

Rat Cd2ap/CD2-associated protein ELISA Kit

Catalog No: #EK19965



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Cd2ap/CD2-associated protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CD2-associated protein
Target Name	CD2-associated protein
Accession No.	F1LRS8
Uniprot	F1LRS8

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.172ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cd2ap

Note: This product is for in vitro research use only