



Mouse Clpp/ATP-dependent Clp protease proteolytic subunit, mitochondrial ELISA Kit

Catalog No: #EK19967

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Clpp/ATP-dependent Clp protease proteolytic subunit, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse ATP-dependent Clp protease proteolytic subunit, mitochondrial
Target Name	ATP-dependent Clp protease proteolytic subunit, mitochondrial
Other Names	Endopeptidase Clp
Accession No.	O88696
Uniprot	O88696
GeneID	53895;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.12ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Clpp

Note: This product is for in vitro research use only