

# Human BCAR3/Breast cancer anti-estrogen resistance protein 3 ELISA Kit

Catalog No: #EK19984

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human BCAR3/Breast cancer anti-estrogen resistance protein 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Breast cancer anti-estrogen resistance protein 3
Target Name	Breast cancer anti-estrogen resistance protein 3
Other Names	Novel SH2-containing protein 2,SH2 domain-containing protein 3B,NSP2,SH2D3B,UNQ271/PRO308
Accession No.	O75815
Uniprot	O75815
GeneID	8412;

## Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.183ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:BCAR3

Note: This product is for in vitro research use only