

Rat Cldn19/Claudin-19 ELISA Kit

Catalog No: #EK19985



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Cldn19/Claudin-19 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Claudin-19
Target Name	Claudin-19
Accession No.	Q5QT56
Uniprot	Q5QT56
GeneID	298487;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.183ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cldn19

Note: This product is for in vitro research use only