

Mouse Adam29/Disintegrin and metalloproteinase domain-containing protein 29 ELISA Kit

Catalog No: #EK19989

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adam29/Disintegrin and metalloproteinase domain-containing protein 29 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Disintegrin and metalloproteinase domain-containing protein 29
Target Name	Disintegrin and metalloproteinase domain-containing protein 29
Other Names	ADAM 29
Accession No.	Q811Q4
Uniprot	Q811Q4

Application Details

Assay Type: Sandwich Detect Range: 0.312-20 ng/mL Sensitivity: 0.151 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: Adam29

Note: This product is for in vitro research use only