

Mouse Cdk12/Cyclin-dependent kinase 12 ELISA Kit

Catalog No: #EK19999



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cdk12/Cyclin-dependent kinase 12 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase 12
Target Name	Cyclin-dependent kinase 12
Other Names	Cdc2-related kinase, arginine/serine-rich,CrkRS,Cell division cycle 2-related protein kinase 7,CDC2-related protein kinase 7,Cell division protein kinase 12,Crk7,Crkrs,Kiaa0904
Accession No.	Q14AX6
Uniprot	Q14AX6
GeneID	69131;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cdk12

Note: This product is for in vitro research use only