

Mouse Angptl7/Angiopoietin-related protein 7 ELISA Kit

Catalog No: #EK20023

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Angptl7/Angiopoietin-related protein 7 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Angiopoietin-related protein 7
Target Name	Angiopoietin-related protein 7
Other Names	Angiopoietin-like protein 7
Accession No.	Q8R1Q3
Uniprot	Q8R1Q3
GeneID	654812;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/ml Sensitivity:0.096ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Angptl7

Note: This product is for in vitro research use only