

## FA58A Antibody

Catalog No: #44720

Package Size: #44720-1 50ul #44720-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FA58A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FA58Antibody detects endogenous levels of total?FA58A
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FA58A
Target Name	FA58A
Other Names	Cyclin-related protein FAM58A,CDK10-activating cyclin,Cyclin-M
Accession No.	Swiss-Prot#: Q8N1B3NCBI Gene ID:92002
Uniprot	Q8N1B3
GeneID	92002;
Calculated MW	29 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



---

Note: This product is for in vitro research use only