

MYOCD Antibody

Catalog No: #44757

Package Size: #44757-1 50ul #44757-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYOCD Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MYOCD antibody detects endogenous levels of total MYOCD
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MYOCD
Target Name	MYOCD
Other Names	MYCD antibody MYCD_HUMAN antibody Myocardin antibody Myocd antibody
Accession No.	Swiss-Prot#: Q8IZQ8 NCBI Gene ID:93649
Uniprot	Q8IZQ8
GeneID	93649;
Calculated MW	100 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western Blot analysis of Hela lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only