

EBP2 Antibody

Catalog No: #44770

Package Size: #44770-1 50ul #44770-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EBP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	EBP2antibody detects endogenous levels of totalEBP2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human EBP2
Target Name	EBP2
Other Names	EBNA1-binding protein 2 antibody ebna1bp2 antibody EBP2 antibody EBP2_HUMAN antibody NOBP antibody nuclear FGF3 binding protein antibody Nucleolar protein p40 antibody P40 antibody Probable rRNA processing protein EBP2 antibody Probable rRNA-processing protein ebp2 antibody
Accession No.	Swiss-Prot#: Q99848NCBI Gene ID:10969
Uniprot	Q99848
GeneID	10969;
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting:□ 1:500-1:2000

Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only