

HSPB3 Antibody

Catalog No: #44808

Package Size: #44808-1 50ul #44808-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | HSPB3 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Immunogen affinity purified |
| Applications | WB |
| Species Reactivity | Hu |
| Specificity | HSPB3antibody detects endogenous levels of total?HSPB3 |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human HSPB3 |
| Target Name | HSPB3 |
| Other Names | 2310035K17Rik antibody A1844863 antibody DHMN2C antibody Heat shock 17 kDa protein antibody Heat shock 27kD protein family, member 3; antibody Heat shock 27kDa protein 3 antibody Heat shock protein 3 antibody Heat shock protein beta 3 antibody HMN2C antibody HSP 17 antibody HSPL27 antibody Protein 3 |
| Accession No. | Swiss-Prot#: Q12988NCBI Gene ID:8988 |
| Uniprot | Q12988 |
| GeneID | 8988; |
| Calculated MW | 17 kDa |
| Concentration | 1.0mg/mL |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Storage | Store at -20°C |

Application Details

Western blotting:□ 1:500-1:2000

Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only