

IGFBPL1 Antibody

Catalog No: #44811

Package Size: #44811-1 50ul #44811-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IGFBPL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	IGFBPL1 antibody detects endogenous levels of total IGFBPL1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IGFBPL1
Target Name	IGFBPL1
Other Names	bA113O24.1 antibody IGFBP-related protein 10 antibody IGFBP-RP4 antibody IGFBPL1 antibody IGFBPRP4 antibody insulin-like growth factor binding protein related protein 4 antibody insulin-like growth factor-binding protein-like 1 antibody insulin-like growth factor-binding-related protein 4 antibody
Accession No.	Swiss-Prot#: Q8WX77NCBI Gene ID:347252
Uniprot	Q8WX77
GeneID	347252;
Calculated MW	29 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only