

MYO1B antibody

Catalog No: #23043

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MYO1B antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 916 and 1125 of MYO1B
Target Name	MYO1B
Accession No.	Swiss-Prot:O43795Gene ID:4430
Uniprot	O43795
GeneID	4430;
Concentration	0.9mg/ml
Formulation	Supplied in 1XPBS, 1%BSA, 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 132kd

Images



Sample (30 ug of whole cell lysate)
A: 293T
7.5% SDS PAGE
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only