

PLCB4 Antibody

Catalog No: #44841



Package Size: #44841-1 50ul #44841-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PLCB4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PLCB4antibody detects endogenous levels of total?PLCB4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PLCB4
Target Name	PLCB4
Other Names	PLCB4
Accession No.	Swiss-Prot#: Q15147NCBI Gene ID:5332
Uniprot	Q15147
GeneID	5332;
Calculated MW	134 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting:□ 1:500-1:2000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only