

platelet factor 4 variant 1 antibody

Catalog No: #23051



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	platelet factor 4 variant 1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 40 and 104 of platelet factor 4 variant 1
Target Name	platelet factor 4 variant 1
Other Names	PF4A; CXCL4L1; CXCL4V1; PF4-ALT; SCYB4V1
Accession No.	Swiss-Prot:P10720Gene ID:5197
Uniprot	P10720
GeneID	5197;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 12kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:250
Immunofluorescence: 1:100-1:200


Images




Sample (30 ug of whole cell lysate)
A: A431
15% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Colon ca, using PF4V1 antibody at 1: 250 dilution.



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed HeLa, using PF4V1 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only