

HTR1D Antibody

Catalog No: #44901

Package Size: #44901-1 50ul #44901-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HTR1D Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	HTR1D antibody detects endogenous levels of total HTR1D
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HTR1D
Target Name	HTR1D
Other Names	5 HT 1D alpha antibody 5 HT 1D antibody 5 HT1B antibody 5 HT1D antibody 5 hydroxytryptamine (serotonin) receptor 1B antibody 5 hydroxytryptamine (serotonin) receptor 1D G protein coupled antibody 5 hydroxytryptamine 1B receptor (5 HT 1B) (Serotonin receptor) (5 HT 1D bet
Accession No.	Swiss-Prot#: P28221NCBI Gene ID:3352
Uniprot	P28221
GeneID	3352;
Calculated MW	45 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of HTR1D expression in HEK293 cells
,The lane on the left is treated with the antigen-specific peptide.

Note: This product is for in vitro research use only