

GPR27 Antibody

Catalog No: #44920

Package Size: #44920-1 50ul #44920-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPR27 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GPR27antibody detects endogenous levels of total?GPR27
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR27
Target Name	GPR27
Other Names	G protein coupled receptor 27 antibody GPR 27 antibody GPR27 antibody GPR27_HUMAN antibody Probable G protein coupled receptor GPR27 antibody Probable G-protein coupled receptor 27 antibody SREB 1 antibody SREB1 antibody Super conserved receptor expressed in brai
Accession No.	Swiss-Prot#: Q9NS67NCBI Gene ID:2850
Uniprot	Q9NS67
GeneID	2850;
Calculated MW	40 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Western blotting:□ 1:500-1:2000


Images



Western blot analysis of lysates from HUVEC cells, using GPR27 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using GPR27 antibody.



Immunofluorescence analysis of A549 cells, using GPR27 Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only