

GPR45 Antibody

Catalog No: #44930

Package Size: #44930-1 50ul #44930-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPR45 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	GPR45antibody detects endogenous levels of total?GPR45
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR45
Target Name	GPR45
Other Names	G protein coupled receptor 45 antibody High affinity lysophosphatidic acid receptor antibody Probable G protein coupled receptor 45 antibody PSP24 1 antibody PSP24 alpha antibody PSP24 antibody PSP24(ALPHA) antibody PSP241 antibody PSP24A antibody
Accession No.	Swiss-Prot#: Q9Y5Y3NCBI Gene ID:11250
Uniprot	Q9Y5Y3
GeneID	11250;
Calculated MW	42 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting:□ 1:500-1:2000

Images



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only