

GPR80/GPR99 Antibody

Catalog No: #44941

Package Size: #44941-1 50ul #44941-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GPR80/GPR99 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GPR80/GPR99antibody detects endogenous levels of total?GPR80/GPR99
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR80/GPR99
Target Name	GPR80/GPR99
Other Names	2 oxoglutarate receptor 1 antibody 2-oxoglutarate receptor 1 antibody aKGR antibody Alpha ketoglutarate receptor 1 antibody Alpha-ketoglutarate receptor 1 antibody G Protein Coupled Receptor 80 antibody G Protein Coupled Receptor 99 antibody G-protein coupled recept
Accession No.	Swiss-Prot#: Q96P68NCBI Gene ID:27199
Uniprot	Q96P68
GeneID	27199;
Calculated MW	38 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting:□ 1:500-1:2000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only