

GPR87/GPR95 Antibody

Catalog No: #44946



Package Size: #44946-1 50ul #44946-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR87/GPR95 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	GPR87/GPR95antibody detects endogenous levels of total?GPR87/GPR95
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR87/GPR95
Target Name	GPR87/GPR95
Other Names	FKSG 78 antibody FKSG78 antibody G protein coupled receptor 87 antibody G protein coupled receptor 95 antibody G-protein coupled receptor 87 antibody G-protein coupled receptor 95 antibody GPCR87 antibody GPR 87 antibody GPR 95 antibody GPR87 antibody GP
Accession No.	Swiss-Prot#: Q9BY21NCBI Gene ID:53836
Uniprot	Q9BY21
GeneID	53836;
Calculated MW	44 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Western blotting:□ 1:500-1:2000

Images



Western Blot analysis of K562 cells using GPR87/95 Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:100

Note: This product is for in vitro research use only