

## OR5H6/5H14 Antibody

Catalog No: #44985

Package Size: #44985-1 50ul #44985-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	OR5H6/5H14 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	OR5H6/5H14 antibody detects endogenous levels of total OR5H6/5H14
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR5H6/5H14
Target Name	OR5H6/5H14
Other Names	Olfactory receptor 5H6 antibody Olfactory receptor OR3-11 antibody OR5H6 antibody OR5H6_HUMAN antibody
Accession No.	Swiss-Prot#: A6NHG9/Q8NGV6NCBI Gene ID:403273/79295
Uniprot	A6NHG9
GenelD	403273;
Calculated MW	37 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western blot analysis of OR5H6/5H14 expression in COLO205 cell lysate

,The lane on the left is treated with the antigen-specific peptide.

---

Note: This product is for in vitro research use only