

## PTGFRN antibody

Catalog No: #23064

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                |
|-----------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name          | PTGFRN antibody                                                                                                |
| Host Species          | Rabbit                                                                                                         |
| Clonality             | Polyclonal                                                                                                     |
| Purification          | Purified by antigen-affinity chromatography.                                                                   |
| Applications          | WB IHC IF                                                                                                      |
| Species Reactivity    | Hu                                                                                                             |
| Immunogen Type        | Peptide                                                                                                        |
| Immunogen Description | Synthetic peptide contain a sequence corresponding to a region within amino acids 815 and 879 of Human PTGFRN  |
| Target Name           | PTGFRN                                                                                                         |
| Other Names           | FPRP; CD315; EWI-F; CD9P-1; SMAP-6; FLJ11001; KIAA1436                                                         |
| Accession No.         | Swiss-Prot:Q9P2B2Gene ID:5738                                                                                  |
| Uniprot               | Q9P2B2                                                                                                         |
| GeneID                | 5738;                                                                                                          |
| Concentration         | 0.7mg/ml                                                                                                       |
| Formulation           | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative. |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                      |

## Application Details


Predicted MW: 99kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Immunofluorescence: 1:100-1:200

## Images




Sample(30 ug of whole cell lysate)


A: H1299

7.5% SDS PAGE


Primary antibody diluted at 1: 500



Immunohistochemical analysis of paraffin-embedded normal oral, using PTGFRN antibody at 1: 100 dilution.



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed A431, using PTGFRN antibody at 1: 500 dilution.

Note: This product is for in vitro research use only