

## FATP1 Antibody

Catalog No: #45112

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FATP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FATP1 Antibody detects endogenous levels of total FATP1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human FATP1
Target Name	FATP1
Other Names	SLC27A1, ACSVL5, FATP, FATP1, FATP-1, Fatty acid transport protein 1
Accession No.	Swiss-Prot#: Q6PCB7NCBI Gene ID:376497
Uniprot	Q6PCB7
GeneID	376497;
Calculated MW	65 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4°over night

Note: This product is for in vitro research use only