

# GCS- alpha-1 Antibody

Catalog No: #45346



Package Size: #45346-1 50ul #45346-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GCS- alpha-1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GCS- alpha-1 Antibody detects endogenous levels of total GCS- alpha-1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GCS- alpha-1
Target Name	GCS- alpha-1
Other Names	GUCY1A3, GC-S-alpha-1, GUC1A3, GUCA3, GCS-alpha-3, GUCY1A1, GC-SA3, GCS-alpha-1, GUCSA3
Accession No.	Swiss-Prot#: Q02108NCBI Gene ID:2982
Uniprot	Q02108
GeneID	2982;
Calculated MW	77 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


## Application Details

Western blotting: 1:1000-3000


## Images



Western blot analysis of lysate from A549 cells, using GCS-α-1 antibody.



Western Blot analysis of various cells using GCS- $\alpha$ -1 Polyclonal Antibody



Immunohistochemistry analysis of GCS- $\alpha$ -1 antibody in paraffin-embedded human lung carcinoma tissue.

Note: This product is for in vitro research use only